

## Animals

Female B<sub>6</sub>C<sub>3</sub>F<sub>1</sub> mice were used in these studies as they are the National Toxicology Program preferred strain for evaluating general toxicity. All mice were purchased from Jackson Laboratory (Bar Harbor, ME) at 7-8 weeks of age. Upon arrival, the animals were allowed to acclimate for a minimum of 5 days. All animals were randomly assigned to treatment groups, weighed, and individually identified via tail marking using a permanent marker. Dose groups were identified by cage cards. Both the dosing group as well as the animal numbers were identified on each cage. The animals were housed 5 mice/cage in ventilated plastic shoe box cages with hardwood chip bedding, modified NIH-31 6% irradiated rodent diet (Harlan Teklad – item #7913) and sterile tap water from water bottles *ad libitum*. The temperature in the animal facility was maintained between 65 and 78 °F and the relative humidity between 30 and 70%; a light/dark cycle was maintained at 12-hr intervals. All animal experiments were performed in an AAALAC International accredited National Institute for Occupational Safety and Health (NIOSH) animal facility in accordance with an animal protocol approved by the Institutional Animal Care and Use Committee.

## Test articles and chemicals

Acetone [CAS #67-64-1] and perfluoropentanoic acid (97%, PFPeA) [CAS #2706-90-3] were purchased from Sigma-Aldrich. Perfluoroheptanoic acid (99%; PFHpA) [CAS# 375-85-9], perfluorohexanoic acid (98%; PFHxA) [CAS# 307-24-4], and perfluorooctanoic acid (95%; PFOA) [CAS# 335-67-1] were purchased from Synquest Laboratories. PFAS concentrations were selected based on an initial 7-day range finding study and previous studies conducted with the short-chain PFAS PFBA (Weatherly et al. 2021) and the long-chain PFAS PFOA (Shane et al. 2020).

## PFAS exposures

For all studies, B<sub>6</sub>C<sub>3</sub>F<sub>1</sub> mice (5/group) were topically treated on the dorsal surface of each ear (25 µl/ear) with vehicle (acetone), PFHpA concentrations (1.25, 2.5, 5% w/v), or PFPeA concentrations (1.25, 2.5, 5% v/v) once a day for 28 consecutive days. Animals were dosed with PFHxA concentrations 2.5, 5, and 10% (v/v) for days 1-4 that was then reduced to 1.25, 2.5, and 5% for days 5-28 due to dermal irritation at the application site. Body weights were measured daily before exposure to ensure no overt toxicity was occurring due to PFAS exposure. Animals were euthanized by CO<sub>2</sub> asphyxiation approximately 24 hrs after the last exposure.

## Tissue processing

Animals were euthanized by CO<sub>2</sub> inhalation, weighed, and examined for gross pathology. The liver, spleen, kidneys, and thymus were removed, cleaned of connective tissue and weighed. Left and right auricular draining lymph nodes (dLNs) and spleen (1/2) were collected in 4 mL RPMI. Spleen (1/2) and dLN (2 nodes/animal) cell suspensions were prepared by mechanical disruption

of tissues between frosted microscope slides in RPMI and cells were counted after RBC lysis using a Z2 Coulter Particle Count and Size Analyzer (Beckman Coulter). One ear pinna was collected and placed in 4 ml of RPMI for immune phenotyping and half of one ear was collected in 0.5ml of RNAlater for subsequent gene expression analysis (see below). Ear cell suspensions were prepared by splitting ear pinna into ventral and dorsal halves, followed by an enzymatic digestion for 90 min at 37 °C with 0.25 mg/ml Liberase-TL Research grade (Roche) in RPMI with 100 µg/ml DNase I (Sigma-Aldrich). Digestion was stopped by the addition of 3ml of RPMI + 10% fetal bovine sera (FBS), the ear pinna + media were transferred to gentleMACS C Tubes (Miltyeni Biotec), then cells were mechanically disrupted on a gentleMACS™ Dissociator (Miltyeni Biotec). Following disruption, cells were passed through a 70 µm cell strainer to make a single cell suspension, washed with RPMI 10% FBS, then live cells were counted on a Cellometer using AO/PI (Nexcelom) in order to quantify cells. A small lobe of the liver was collected in 0.5 ml of RNAlater for subsequent gene expression (see below). The remainder of the liver, spleen (1/2), ear (1/2), and one kidney was collected in 10% formalin for histopathology analysis.

### **Serum chemistries**

Blood samples were collected via cardiac puncture, transferred to serum separation tubes, and separated by centrifugation. The serum was frozen at -20 °C for subsequent serum chemistry analysis. Selected serum chemistries were evaluated using a Catalyst DX Chemistry Analyzer (IDEXX Laboratories, Inc.; Westbrook, ME). Endpoints analyzed included: alkaline phosphates (ALKP), urea nitrogen (BUN), glucose (GLU), and cholesterol (CHOL).

### **Analytical PFAS detection**

Serum and urine samples were collected and analyzed for individual PFAS (PFPeA, PFHxA, PFHpA) by Vista Analytical Laboratory, following Vista's standard operating procedures (Vista Analytical Laboratory, El Dorado Hills, CA) of solid phase extraction and liquid chromatography/tandem mass spectrometry (LC/MS/MS) as described in (Shoemaker et al. 2008). Serum was collected from each animal and urine was pooled following euthanasia for each group of mice. The results for PFPeA, PFHxA, and PFHpA include the linear isomer only. The Initial Calibration and Continuing Calibration Verifications met the acceptance criteria as described in (Shoemaker et al. 2008). No analytes were detected in the Method Blank above Reporting Limit. The labeled standard recoveries for all quality controls and samples were within the acceptance criteria as described in (Shoemaker et al. 2008).

### **Flow Cytometry**

For staining, single cell suspensions were resuspended in staining buffer containing α-mouse CD16/32 antibody (Fc Block) (BD Biosciences) then incubated with a cocktail of fluorochrome-conjugated antibodies specific for mouse cell surface antigens. For dLN and spleen cells: B220-V500 (RA3-6B2), CD11b-PerCP-Cy5.5 (M1/70), CD8-PE-CF594 (53-6.7), Siglec-F-PE (E50-

2440), Ly6G-FITC (1A8) (BD Biosciences), CD11c-eFluor 450 (N418), CD86-APC (GL1) (eBioscience), CD45-Superbright 780 (30-F11), F4/80-PE-Cy7 (BM8), MHCII-APC-eF780 (M5/114.15.2) (Invitrogen), CD4-BV711 (RM4-5), Ly6C-AF700 (HK1.4) (BioLegend). For ear cells: CD4-BV711 (RM4-5), NKp46-BV605 (29A1.4) (BioLegend), CD8-PE-CF594 (53-6.7), CD3-V500 (500A2), Ly6G-FITC (1A8), CD11b-PerCP-Cy5.5 (M1/70), CD11c-AF700 (HL3), Siglec-F-PE (E50-2440) (BD Biosciences), F4/80-PE-Cy7 (BM8), MHCII-APC-eF780 (M5/114.15.2), CD45-Superbright 780 (30-F11), FcεRI-APC (MAR-1) (Invitrogen), CD117-eFluor 450 (2B8) (eBioscience). Cells were then washed, fixed in Cytofix buffer (BD Biosciences), resuspended in staining buffer, and a minimum of 100,000 CD45<sup>+</sup> events were collected on a LSR II flow cytometer (BD Biosciences), compensation controls were prepared with eBioscience UltraComp eBeads. Analysis was performed using FlowJo v10 software (TreeStar Inc., Ashland, OR). All events were gated on single cells using FSC and SSC parameters prior to subsequent gating. Cellular populations were defined using the gating strategies outlined below; Fluorescence minus ones (FMOs) were used as gating controls.

### Gating strategies for immunophenotyping

Spleen and dLN cells	
Cell Type	Gating Strategy
CD4 T cells	CD45 <sup>+</sup> , lymphocyte gate (FSC/SSC <sup>low</sup> ), CD8 <sup>-</sup> , B220 <sup>-</sup> , NKp46 <sup>-</sup> , CD4 <sup>+</sup>
CD8 T cells	CD45 <sup>+</sup> , lymphocyte gate (FSC/SSC <sup>low</sup> ), CD4 <sup>-</sup> , B220 <sup>-</sup> , NKp46 <sup>-</sup> , CD8 <sup>+</sup>
B cells	CD45 <sup>+</sup> , lymphocyte gate (FSC/SSC <sup>low</sup> ), CD4 <sup>-</sup> , CD8 <sup>-</sup> , NKp46 <sup>-</sup> , B220 <sup>+</sup> (MHCII/CD86 MFI increases with activation)
NK cells	CD45 <sup>+</sup> , lymphocyte gate (FSC/SSC <sup>low</sup> ), CD4 <sup>-</sup> , CD8 <sup>-</sup> , B220 <sup>-</sup> , NKp46 <sup>+</sup>
Eosinophils	CD45 <sup>+</sup> , SSC <sup>high</sup> , CD11b <sup>+</sup> , Ly6G <sup>low</sup> , Siglec-F <sup>+</sup>
Neutrophils	CD45 <sup>+</sup> , SSC <sup>high</sup> , CD11b <sup>+</sup> , Siglec-F <sup>-</sup> , Ly6G <sup>high</sup>
DCs	CD45 <sup>+</sup> , SSC <sup>int</sup> , MHCII <sup>+</sup> , CD11c <sup>+</sup> , F4/80 <sup>-</sup> (MHCII/CD86 MFI increases with activation)
CD11b <sup>+</sup> monocytes/macrophages	CD45 <sup>+</sup> , SSC <sup>int</sup> , CD11b <sup>+</sup> , SiglecF <sup>-</sup> , Ly6G <sup>-</sup> , Lyc6 <sup>+/-</sup> , F4/80 <sup>+/-</sup>
Skin (ear) cells	
Cell Type	Gating Strategy
CD4 T cells	CD45 <sup>+</sup> , Lymphocyte gate (FSC/SSC <sup>low</sup> ), CD3 <sup>+</sup> , CD4 <sup>+</sup> , CD8 <sup>-</sup>
CD8 T cells	CD45 <sup>+</sup> , Lymphocyte gate (FSC/SSC <sup>low</sup> ), CD3 <sup>+</sup> , CD4 <sup>-</sup> , CD8 <sup>+</sup>
NK cells	CD45 <sup>+</sup> , Lymphocyte gate (FSC/SSC <sup>low</sup> ), CD3 <sup>-</sup> , NKp46 <sup>+</sup>
Eosinophils	CD45 <sup>+</sup> , SSC <sup>high</sup> , CD11b <sup>+</sup> , Ly6G <sup>low</sup> , Siglec-F <sup>+</sup>
Neutrophils	CD45 <sup>+</sup> , SSC <sup>high</sup> , CD11b <sup>+</sup> , Siglec-F <sup>-</sup> , Ly6G <sup>high</sup>
DCs	CD45 <sup>+</sup> , SSC <sup>int</sup> , MHCII <sup>+</sup> , CD11c <sup>+</sup> , F4/80 <sup>-</sup> , CD11b <sup>+/-</sup>

## Gene expression

Ear (half an ear pinna/mouse) and liver (caudate) were homogenized on a TissueLyser II in Buffer RLT (Qiagen). Total RNA was isolated using Qiagen's RNeasy mini spin column kits with DNase treatment on a QIAcube automated RNA isolation machine. RNA concentrations and purity were analyzed on a NanoDrop spectrophotometer (Thermo Fisher Scientific). The cDNA (1-2 µg) was prepared on an Eppendorf Mastercycler using Applied Biosystems' High Capacity Reverse Transcription kit. The cDNA was used as a template for real-time PCR reactions containing TaqMan PCR Master Mix with gene-specific primers (Applied Biosystems) on a 7500 Real-Time PCR System. Relative fold gene expression changes ( $2^{-\Delta\Delta CT}$ ) were determined compared to vehicle controls and normalized for expression of reference gene beta-actin (Taqman).

## Histology

Tissues were collected in 10% formalin for histopathology. Each tissue sample was embedded in paraffin, sectioned at 5 microns, stained with hematoxylin and eosin (H&E) and evaluated by a veterinary pathologist at StageBio (Mason, Ohio) using The Society of Toxicologic Pathology Guideline (Crissman et al. 2004). Provantis™ pathology software v10.2.3.1 was utilized for data capture and table generation. Histopathology grades were assigned as grade 1 (minimal), grade 2 (mild), grade 3 (moderate), grade 4 (marked), or grade 5 (severe) based on an increasing extent of change.

## Histopathology grading criteria

### **Dosing Site Scoring**

#### Epidermal hyperplasia

- |             |                             |
|-------------|-----------------------------|
| 0- Normal   | 1-2 layers of keratinocytes |
| 1- Minimal  | 3-4 layers of keratinocytes |
| 2- Mild     | 5-6 layers of keratinocytes |
| 3- Moderate | >6 layers of keratinocytes  |

#### Dermal Fibrosis

- |              |                                      |
|--------------|--------------------------------------|
| 0- Normal    |                                      |
| 1 – Minimal  | Increased thickness, but <2X control |
| 2- Mild      | 2-3X control thickness               |
| 3 – Moderate | >3X control thickness                |

#### Epidermal Hyperkeratosis

- |                   |  |
|-------------------|--|
| 0- Absent         |  |
| 1- Minimal –      | <2X control of compacted keratin         |
| 2- Mild - approx. | 2X control amounts of compacted keratin  |
| 3- Moderate -     | >2X control amounts of compacted keratin |

#### Erosion/ulcer (associated with inflammation only)

- |              |           |
|--------------|-----------|
| 0- Absent    |           |
| 1- Minimal – | One focus |
| 2- Mild -    | Two foci. |

#### Dermal Inflammation

- |              |  |
|--------------|--|
| 0- Normal -  | No increase in leukocytes compared to controls   |
| 1- Minimal - | Diffuse minimal increase in leukocyte density scattered in dermis or small focus associated with erosion/ulcer |
| 2- Mild -    | Above, with small accumulations of inflammatory cells in the dermis  |

#### Epidermal necrosis (coagulative)

- |               |   |
|---------------|---|
| 0- Absent     |   |
| 1- Minimal -  | Superficial loss of differential staining affecting <10% of dosed surface of ear.   |
| 2- Mild -     | Loss of differential staining affecting <10% of dosed side of ear, full thickness in at least one area or superficial loss of differential staining for >10% of surface area. |
| 3- Moderate - | Loss of differential staining affecting ~10-40% of surface area with full thickness in at least one area.   |
| 4- Marked -   | Loss of differential staining affecting ~40-90% of surface area with full thickness in at least one area.   |

- Crissman JW, Goodman DG, Hildebrandt PK, Maronpot RR, Prater DA, Riley JH, Seaman WJ, Thake DC. 2004. Best practices guideline: Toxicologic histopathology. *Toxicologic pathology*. 32(1):126-131.
- Shane HL, Baur R, Lukomska E, Weatherly L, Anderson SE. 2020. Immunotoxicity and allergenic potential induced by topical application of perfluorooctanoic acid (pfoa) in a murine model. *Food and chemical toxicology : an international journal published for the British Industrial Biological Research Association*. 136:111114.
- Shoemaker JA, Grimmett P, Boutin B. 2008. Determination of selected perfluorinated alkyl acids in drinking water by solid phase extraction and liquid chromatography/tandem mass spectrometry (lc/ms/ms).  
[https://cfpub.epa.gov/si/si\\_public\\_record\\_report.cfm?Lab=NERL&dirEntryId=198984&simpleSearch=1&searchAll=EPA%2F600%2FR-08%2F092+](https://cfpub.epa.gov/si/si_public_record_report.cfm?Lab=NERL&dirEntryId=198984&simpleSearch=1&searchAll=EPA%2F600%2FR-08%2F092+).
- Weatherly LM, Shane HL, Lukomska E, Baur R, Anderson SE. 2021. Systemic toxicity induced by topical application of heptafluorobutyric acid (pfba) in a murine model. *Food and chemical toxicology : an international journal published for the British Industrial Biological Research Association*. 156:112528.